# Enter a New Era in Large Bore Closure

## MANTA Vascular Closure Device





MANTA

# A new era of simplicity and confidence in closure is here

The MANTA Device is the first commercially available biomechanical vascular closure device designed specifically for large bore femoral arterial access site closure.<sup>1</sup>

# Simple deployment

Addresses the challenges of large bore closure with a single easy-to-use device.<sup>2a</sup>

## **Rapid** hemostasis

Reduces time to hemostasis without pre-closure, utilizing the coagulation-inducing properties of collagen for rapid hemostasis to promote vessel healing.<sup>2b,3-5</sup> مر المراب ال

## **Reliable** closure

Delivers reproducible results and helps inspire confidence in achieving successful closure.<sup>2c</sup>

> **Over-the-wire design** aids device placement throughout deployment

#### 1. Data on file at Teleflex.

- 2. Data on file at Teleflex. The SAFE MANTA IDE Clinical Trial.
- a. A single MANTA Vascular Closure Device was deployed in 99.6% of subjects in IDE trial.
- b. The MANTA Device demonstrated a time to hemostasis (TTH) of 24 seconds median time (65 seconds mean time) from deployment to hemostasis, which is lower than published
- c. 97.7% Technical Success, defined as percutaneous vascular closure obtained with the MANTA Device without the use of unplanned endovascular or surgical intervention.
- Study sponsored by Teleflex Incorporated or its affiliates.
- J Vasc Surg. 2014 May;59(5):1081-1193
- 4. Farndale RW, et al. The role of collagen in thrombosis and hemostasis. J Thromb Haemost. 2004 Apr;2(4):564-573.
- 5. Nuyttens BP, et al. Platelet adhesion to collagen. Thromb Res. 2011;127(2): S26-S29.
- a. Percutaneous vascular closure obtained with the MANTA Device without the use of unplanned endovascular or surgical intervention.
- b. The MANTA Device demonstrated a time to hemostasis (TTH) of 24 seconds median time (65 seconds mean time) from deployment to hemostasis.
- c. Major complications defined as composite of i) vascular injury requiring surgical repair/stent-graft; ii) bleeding requiring transfusion; iii) lower extremity ischemia requiring surgical repair/additional percutaneous intervention; iv) nerve injury (permanent or requiring surgical repair); and v) infection requiring IV antibiotics and/or extended hospitalization Study sponsored by Teleflex Incorporated or its affiliates.
- 8. Lauten A, et al. Percutaneous left-ventricular support with the Impella 2.5°-assist device in acute cardiogenic
- shock: results of the Impella-EUROSHOCK-registry. Circ Heart Fail 2013 Jan;6(1):23-30



### **Engineered for versatility**

Available in 14 Fr. and 18 Fr., a single MANTA Device effectively closes femoral arterial access sites following the use of large bore sheaths ranging from 12 Fr. to 25 Fr. O.D.

### Ordering Information

#### The Manta Vascular Closure Device

The 14F MANTA is indicated for closure of femoral arterial access sites following the use of 10-14F devices or sheaths (maximum OD/profile of 18F), and the 18F MANTA device is indicated for closure of femoral arterial access sites following the use of 15-18F devices or sheaths (maximum OD/profile of 25F).

MODEL	DESCRIPTION	SIZE
2156NE	14 Fr. MANTA Vascular Closure Device	14Fr
2115NE	18 Fr. MANTA Vascular Closure Device	18Fr

Packaged in quantities of 5 unit per box.

# How it works

The MANTA Device facilitates biomechanical closure without pre-closure.



1. Insert the MANTA Device







3. Withdraw and seal



### **Front view**

secures the sliding suture knot without tamping and is a helpful landmark for future interventions

#### Sliding suture knot

provides initial compaction

## **Resorbable collagen**

sandwich the access site

# Clinically proven

Tale

13



patients

263

The SAFE MANTA IDE Clinical Trial- the largest U.S. prospective multi-center study of a purpose-designed large bore femoral arterial access site closure device to date -demonstrated the safety and effectiveness of the MANTA Device.<sup>6</sup>

MANTA



Technical success rate<sup>6a</sup>



Median time from deployment to hemostasis (65 seconds mean time)<sup>6b</sup>

 $5.3^{0/0}$ 

#### Major Complications rate

Breakdown by most serious event or intervention: Major Bleeding 2.3% (6/263) Covered Stent 1.5% (4/263) Balloon 0.8% (2/263) Surgery 0.8% (2/263)



intervention: Covered Stent 1.5% (4/263) Major Bleeding 1.1% (3/263) Balloon 0.8% (2/263) Surgery 0.8% (2/263) **INDICATIONS:** The 14F MANTA is indicated for closure of femoral arterial access sites following the use of 10-14F devices or sheaths (maximum OD/ profile or 18F), and the 18F MANTA device is indicated for closure of femoral arterial access sites following the use of 15-18F devices or sheaths (maximum OD/profile of 25F).

**CONTRAINDICATIONS:** 1) Severe calcification of the access vessel; 2) Severe peripheral artery disease; 3) Puncture in the origin of the profundal femoral artery, above the inguinal ligament, or above the most inferior border of the epigastric artery (IEA); 4) Sheath insertion in vessel other than the femoral artery; 5) Marked tortuosity of the femoral or iliac artery; 6) Marked obesity or cachexia (BMI >40 or <20); 7) Blood pressure >180 mmHg; 8) Patients who cannot be anti-coagulated for the procedure.

WARNINGS: Do not use: 1) if bacterial contamination of procedure sheath or surrounding tissues may have occurred; 2) if the procedure sheath has been placed through the superficial femoral artery and into the profunda femoris artery; 3) if the MANTA Device delivery system becomes kinked; 4) with a contralateral balloon inflated in the femoral or iliac artery during MANTA device use; 5) if there has been a femoral artery puncture in same vessel within prior 30 days, recent femoral artery puncture in same groin that has not healed appropriately, and/or recent (<30 days) vascular closure device placement in same femoral artery; 6) if puncture site is at or distal to bifurcation of superficial femoral and profunda femoris artery; 7) if the puncture site is proximal to inguinal ligament or above most inferior border of epigastric artery (IEA).

**POTENTIAL ADVERSE EVENTS:** 1) Failed hemostasis requiring manual or mechanical compression, application of balloon pressure from a secondary access site, placement of a covered stent or surgical repair. 2) Local trauma to the femoral or iliac artery wall, such as dissection. 3) Retroperitoneal bleeding as a result of access above the inguinal ligament or the most inferior border of the epigastric artery (IEA). 4) Perforation of ileofemoral arteries, causing bleeding/hemorrhage. 5) Accidental positioning of some or all of the collagen plug within the femoral artery, leading to ischemia or stenosis. 6) Thrombosis formation or embolism. 7) Nerve damage or neuropathy. 8) Other access site complications leading to bleeding, hematoma, pseudoaneurysm, etc., possibly requiring blood transfusion, surgical repair, and/or endovascular intervention.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose-driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow, Deknatel, Hudson RCI, LMA, Pilling, Rüsch, UroLift and Weck – trusted brands united by a common sense of purpose.

#### **Corporate Office**

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

#### **Regional Offices**

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Teleflex, the Teleflex logo, Arrow, Deknatel, Hudson RCI, LMA, MANTA, Pilling, Rüsch, UroLift and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks and registered trademarks are property of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 02/2019.

© 2019 Teleflex Incorporated. All rights reserved.

MCI-2019-0235-EN REV 1 06 19 PDF

